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[www.finescience.de](http://www.finescience.de)  
[europe@finescience.de](mailto:europe@finescience.de)

**DEUTSCH** (ab S. 1)

**ENGLISH** (from p. 6)

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**ITALIANO** (da p. 16)

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### Instruction Manual

#### **CAUTION!**

The TR-200 is not compatible with older Tissue Slice Chamber No. 21000-02. Please contact FST if you wish to replace your old TR-100.

Thank you for purchasing the TR-200 Temperature Controller (FST 21052-00). The following items are enclosed:

- One TR-200 Temperature Controller
- One Power Supply Cable with two DB-25 connectors

If you also purchased a power supply unit from FST (FST 21051-00), you will find the additional item enclosed:

- One 24 VDC ( $\pm 1$ V), 4 Amp Power Supply

### **Connecting the TR-200 to the Tissue Slice Recording Chamber**

The power requirements for the TR-200 are 24 VDC, 5 A. Two diagrams are provided with this documentation, showing the front and back panels of the TR-200 (see page 5). The various connections are identified and numbered.

#### **1. Connection of the FST Power Supply (No 21051-00) to the TR-200**

With the power switch in the off position, connect one end of the Power Supply Cord to the Power Supply Unit and the other end to the buildings main outlet. Assure that the correct main voltage is showing on the red slide (either 115V or 230V). The Power Supply Input Connector is located on the rear panel of the Temperature Controller (figure 1; number 1). The female DB-25 connector attaches to the corresponding Power Supply Connector and the male DB-25 connector attaches to the TR-200. Do not attempt to loosen the DB-25 connector screws on the bodies of the Power Supply Unit or the TR -200, as they bolt the connectors to their housing.

#### **2. Connection of the TR-200 to the Tissue Slice Recording Chamber**

Connect the leads of the orange Heating Blanket in the water bath of the Tissue Slice Recording Chamber to the white connector of the Heating Blanket Cable. Next, insert the banana plugs of the Heating Blanket Cable into the power outlet on the rear panel of the TR-200 (figure 1, number 3). Note, that there is no polarity to the Heating Element.

There are three thermistors in the Tissue Slice Recording Chamber. Since the temperature in each recording well will be nearly identical, we recommend that you connect one set of the thermistor leads to the bath thermistor (located at the top of the recording head) and the other set of thermistor leads to one of the recording well thermistors. This will allow you to control the temperature of the large water bath, and read the temperature of the recording well, by setting the Thermistor Selection Switch (figure 2, number 6) to position "A" (suggested to be the water bath thermistor) or position "B". The TR-200 is calibrated to use the thermistors that are incorporated into the Tissue Slice Recording Chamber which have a resistance of 25ohms at 25 degrees Celsius.

### **Operation of the TR-200 and the Tissue Slice Recording Chamber**

*Note: Please refer to the Instruction Manual of the Tissue Slice Recording Chamber that was provided with that unit.*

Once you have connected the Power Supply, the TR-200 and the Tissue Slice Recording Chamber, it is now a simple matter to begin the maintenance of the temperature of your tissue slices. First, fill the chamber of the water bath with water (see the instruction booklet for the Tissue Slice Recording Chamber). Make sure that both the Power Supply and the TR-200 are turned on, and the Thermistor Selection Switch is in the “A” position. The TR-200 is turned on by pulling the Temperature Control Dial (figure 2, number 8) until you hear a click. The low initial temperature being measured by the bath thermistor should be displayed on the Temperature Display (Figure 2, number 4) in degrees Celsius.

To set the temperature at which you want your tissue slices to be maintained, push down the Run/Set Switch (figure 2, number 5) and, while holding it down, turn the Temperature Control Dial (figure 2, number 8) clockwise until the desired maintenance temperature is displayed on the Temperature Display. Now release the Run/Set Switch. If the temperature in the Tissue Slice Recording Chamber is lower than the temperature that you have set (which it is likely to be), the red Power Indicator Light (figure 2, number 9) will come on to indicate that power is being supplied to the Heating Blanket. Once the water bath has reached the desired temperature, then the power to the Heating Blanket will stop and the green Power Indicator Light will come on.

**CAUTION:** If you turn the TR-200 off, the pre-set temperature at which you had your tissue slices maintained will have changed. You will have to re-set the desired temperature again for your next experiment. ALSO, once you have set the desired temperature, DO NOT turn the Temperature Control Dial without depressing the Run/Set button, as this will invalidate the previously set temperature.

### **Calibration of the TR-200**

The Temperature Controller is pre-calibrated by the manufacturer prior to shipment. If, although unlikely, it becomes necessary to re-calibrate the TR-200, adjust the two calibration resistors by following procedures 1 and 2 below.

#### **Procedure 1: The Offset Adjustment**

The TR-200 is calibrated by way of two fine screw offset adjustment resistors. They are accessed by removing the top cover of the temperature controller. Remove the four screws at the sides of the unit and you find the OFFSET adjusting resistor to be the left one of the two facing up on the circuit board. The procedure is as follows:

Attach one of the Thermistors that are supplied with the Tissue Slice Recording Chamber into its respective input on the front of the temperature controller. Make sure that the Thermistor Selector Switch is pointing to the correct side. Place the tip of the thermistor into a water bath of known temperature, (using a thermometer accurate to the nearest tenth of a degree Celsius). If the Temperature Display does not correspond to the thermometer, turn the screw of the left resistor until the Temperature Display corresponds to the reading of the thermometer.

## Procedure 2: Gain Adjustment

The second calibration resistor is also accessed by removing the top cover of the Temperature Controller. The GAIN adjustment resistor is the right one of the two fine screw resistors on the main circuit board. To adjust the gain, follow this procedure:

After calibrating the offset adjustment (procedure 1 above) set the water bath to a higher temperature than before. If the Temperature Display does not agree with the thermometer reading a gain adjustment is needed. In order to increase the gain, turn the screw on the right resistor clockwise until both temperatures are equal; to decrease the gain, turn the screw counterclockwise. For better accuracy the procedures should be done with a 10 degree temperature difference. For example, perform procedure 1 at 30°C and procedure 2 at 40°C.

## Fuse Replacement

To open the fuse container, push the fuse container cap (figure 1, number 2) inward and rotate counterclockwise to remove. Pull fuse from socket and replace with a new 5A to 6A, 125V fuse. Align fuse holder and push into socket. Push against spring and rotate cap clockwise.

## Diagrams


Figure 1: Rear Panel of the Temperature Controller


Figure 2: Front Panel of the Temperature Controller

## **Part Descriptions**

1. Power Supply Input Connector
2. Fuse Container
3. Tissue Slice Recording Chamber Power Output
4. Temperature Display
5. Run/Set Switch
6. Thermistor Selection
7. Thermistor Probe Inputs
8. Temperature Control Dial
9. Power Indicator Lights

\* Accuracy is +/- 0.2° C in an environment of 1V/m field strength and +/- 1°C in an environment of 3V/m field strength.

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