


Excellence Analytical Balances

XS Models – Part 1



METTLER TOLEDO

Table of contents

1	Introduction		5
	1.1	Conventions and Symbols Used in These Operating Instructions	6
2	Safety Information		7
	2.1	Definition of Signal Warnings and Symbols	7
	2.2	Product Specific Safety Notes	7
3	Overview XS Analytical Balance		9
4	Setting up the Balance		10
	4.1	Unpacking and Checking the Delivered Items	10
	4.1.1	Unpacking the Balance	10
	4.2	Scope of Delivery	11
	4.3	Selecting a Location and Leveling the Balance	12
	4.3.1	Location	12
	4.3.2	Leveling the Balance	12
	4.4	Assembling the Balance	13
	4.5	Power Supply	15
	4.6	Left/Right Operating of the Glass Draft Shield	16
	4.7	Setting the Reading Angle and Positioning the Terminal	16
	4.7.1	Optimise the Readability of the Terminal	16
	4.7.2	Remove Terminal and Place Close to the Balance	16
	4.8	Transporting the Balance	17
	4.8.1	Transporting Over Short Distances	17
	4.8.2	Transporting Over Long Distances	18
	4.9	Below-the-Balance Weighing	20
	4.10	Installing the ErgoClip	20
	4.11	Installing the Single-use Aluminum Weighing Pan	21
	4.12	Installing the Grid Weighing Pan Cover	21
5	Maintenance		23
	5.1	Cleaning	23
	5.2	Disposal	23
6	Technical Data		24
	6.1	General Data	24
	6.2	Explanatory Notes for the METTLER TOLEDO AC Adapter	24
	6.3	Model-specific Data	25
	6.4	Dimensions	28
	6.5	Interfaces	29
	6.5.1	Specifications of RS232C	29
	6.5.2	Specifications of "Aux" Connection	29
7	Accessories and Spare Parts		30
	7.1	Accessories	30
	7.2	Spare Parts	41

8	Appendix		42	
		8.1	MT-SICS Interface Commands and Functions	42
		8.2	Procedure for Certified Balances	42
9	Index		44	

1 Introduction

Thank you for choosing a METTLER TOLEDO balance.

The balances of the XS line combine a large number of weighing and adjustment possibilities with exceptionally convenient operation.

In this chapter you will be given basic information about your balance. Please read right through this chapter carefully even if you already have experience with METTLER TOLEDO balances. Please pay special attention to the safety warnings!

The different models have different characteristics regarding equipment and performance. Special notes in the text indicate where this makes a difference to operation.

The XS line comprises a range of balances which differ from each other in relation to their weighing range and resolution.

The following features are common to all models of the XS lines:

- Fully automatic adjustment "FACT" using internal weight.
- Built-in applications for normal weighing, statistics, formulation, density, percent weighing, piece counting and LabX Client.
- Integral RS232C interface.
- Slot for second interface (optional).
- Touch-sensitive graphics terminal ("Touch Screen") for easy, convenient operation.

A brief word about standards, guidelines, and methods of quality assurance: The balances comply with usual standards and guidelines. They support standard procedures, specifications, working methods, and reports according to **GLP (Good Laboratory Practice)**. In this connection, records of working procedures and adjustments become very important; for this purpose we recommend you to use a printer from the METTLER TOLEDO range, since these are optimally adapted to your balance. The balances conform to the applicable standards and guidelines and possess a EC declaration of conformity. METTLER TOLEDO is certified as manufacturer according to ISO 9001 and ISO 14001.

The Operating Instructions for the XS balances consist of 3 separate documents, whose contents are listed below.

Part 1, This Document

Contents

- Introduction
- Safety Information
- Setting up the Balance
- Leveling the Balance
- Cleaning and Service
- Technical Data
- Interface commands and MT-SICS functions
- Accessories
- Spare Parts

Part 2, Separate Document

Contents: Terminal, System and Applications

- Basic Principles for Using the Terminal and the Firmware
- System Settings
- Applications
- Firmware (Software) Updates
- Error and Status Messages

- Conversion Table for Weight Units
- Recommended Printer Settings

Part 3, Separate Document

Contents: Adjustments and Tests

- Adjustments
- Tests

Finding More Information

Internet <http://www.mt.com/excellence>

1.1 Conventions and Symbols Used in These Operating Instructions

The following conventions apply to the operating instructions: Part 1, Part 2 and Part 3.

Key designations are indicated by a picture or text in double angular parentheses (e.g. «» or «**On/Off**»).



This symbol indicates press key briefly (less than 1.5 s).



This symbol indicates press and hold key down (longer than 1.5 s).

These symbols indicate an instruction:

- ▶ prerequisites
- 1 steps
- 2 ...
- ⇒ results

2 Safety Information

2.1 Definition of Signal Warnings and Symbols

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the balance, malfunctions and false results.

Signal Words

WARNING	for a hazardous situation with medium risk, possibly resulting in severe injuries or death if not avoided.
CAUTION	for a hazardous situation with low risk, resulting in damaged to the device or the property or in losing of data or minor or medium injuries if not avoided.
Attention	(no symbol) for important information about the product.
Note	(no symbol) for useful information about the product.

Warning Symbols



General hazard



Electrical shock

2.2 Product Specific Safety Notes

Always operate and use your balance only in accordance with the Operating Instructions Part 1, Part 2 and Part 3.

The instructions for setting up your new balance must be strictly observed.

If the instrument is not used according to the manufacturer's Operating Instructions (Part 1, Part 2 and Part 3), protection of the instrument may be impaired.

Intended Use

Your balance is used for weighing. Use the balance exclusively for this purpose. Any other type of use and operation beyond the limits of technical specifications without written consent from Mettler-Toledo AG, is considered as not intended.



It is not permitted to use the instrument in explosive atmosphere of gases, steam, fog, dust and flammable dust (hazardous environments).



CAUTION

Damage of Device

- For use only in dry interior rooms.
 - Do not use sharply pointed objects to operate the keyboard! Although your balance is very ruggedly constructed, it is nevertheless a precision instrument. Treat it with corresponding care.
 - Do not open the balance: It does not contain any parts which can be maintained, repaired, or replaced by the user. If you ever have problems with your balance, contact your METTLER TOLEDO dealer.
 - Use only balance accessories and peripheral devices from METTLER TOLEDO; they are optimally adapted to your balance.
-




CAUTION

Damage of Device

Use only the original universal AC adapter delivered with your balance, and check that the voltage printed on it is the same as your local power supply voltage. Only plug the adapter into a socket which is grounded.

3 Overview XS Analytical Balance



Overview

1	Terminal (details see Operating Instructions – Part 2)	2	Display (Touch-sensitive "Touch Screen")
3	Operating keys	4	Drip tray
5	Handle for variable operation of the draft-shield side doors	6	Grid weighing pan
7	Glass draft shield	8	Type name
9	Handle for operation of the draft-shield top door	10	Guide for top door of draft shield and handle for transport
11	Level indicator	12	Foot screw
13	Slot for second interface (optional)	14	Socket for AC adapter
15	Fastening point for anti-theft device	16	RS232C serial interface
17	Aux 2 (connection for "ErgoSens", hand- or foot-switch)	18	Aux 1 (connection for "ErgoSens", hand- or foot-switch)

4 Setting up the Balance

This chapter explains how to unpack your new balance, and how to set it up and prepare it for operation. When you have carried out the steps described in this chapter, your balance is ready for operation.

4.1 Unpacking and Checking the Delivered Items

4.1.1 Unpacking the Balance

Use the lifting strap to lift the balance out of the packaging carton.


Overview

- 1 Lifting strap
- 2 Top packing cushion
- 3 Operating Instructions and other important documents
- 4 Balance
- 5 Set with draft-shield doors and terminal support
- 6 Set with AC adapter, power supply cable, drip tray, grid weighing pan, grid weighing pan cover, set of single-use pans and Ergo-Clip "Basket" (basket for small weighing objects)
- 7 Terminal
- 8 Bottom packing cushion


Note

The terminal is connected to the balance by a cable!

- 8 Bottom packing cushion



- 1 Unfasten lifting strap (1).
- 2 Remove top packing cushion (2).




- 1 Pull out Operating Instructions (3).
- 2 Remove set with AC adapter etc. (4).
- 3 Remove set with draft-shield doors etc. (5).




- Carefully pull the terminal out of the bottom packing cushion and remove the protective cover.

Note


The terminal is connected to the balance by a cable, so only pull the terminal just far enough out of the packing cushion to remove the protective cover.



- 1 Place the terminal on the front of the balance.
- 2 Hold the balance by the guide or handle, hold the terminal firmly with your other hand, and pull the balance and terminal together out of the bottom packing cushion.



- 1 Place the balance with the terminal in the place where the balance will be used for weighing.
- 2 Remove the cover from the balance.
- 3 Pull the transport protection (9) of the weighing pan support toward the front and off.



Note

Please keep all parts of the packaging. This packaging guarantees best possible protection of your balance for transportation, **see** Transporting the Balance (page 17).

4.2 Scope of Delivery

The standard scope of delivery contains the following items:

- Balance with terminal
 - RS232C interface
 - Slot for second interface (optional)
 - Feedthroughs for below-the-balance weighing and for anti-theft device

- Set with draft-shield doors and terminal support
- Grid weighing pan
- Grid weighing pan cover of chrome-nickel steel (attachment for grid weighing pan)
- Set of single-use aluminum weighing pans (10 pans) for mounting on the grid weighing pan
- Drip tray
- AC adapter with country-specific power cable
- Protective cover for the terminal
- Cleaning brush
- ErgoClip "Basket" (basket for small weighing objects)
- Production certificate
- CE declaration of conformity
- Operating Instructions Part 1 (this document), Part 2 and Part 3
- Instructions for unpacking, re-packing, and setting up

4.3 Selecting a Location and Leveling the Balance


4.3.1 Location

Choose a position which is stable, free from vibration, and as nearly horizontal as possible. The supporting surface must be able to bear the weight of the fully loaded balance safely.

Avoid the following:

- Direct sunlight
- Draft (e.g. from fans or air conditioning)
- Excessive fluctuations in temperature.

Further information can be found in Weighing the Right Way.



Observe ambient conditions. **See** Technical Data (page 24).

4.3.2 Leveling the Balance

The balances have a level indicator and two adjustable leveling feet to compensate for slight irregularities in the surface of the weighing bench.

- The balance is exactly horizontal when the air bubble is in the middle of the level glass.

- 1 Adjust the two leveling feet appropriately until the air bubble comes to rest exactly in the middle of the glass:


Air bubble at "12 o'clock" turn both feet counterclockwise.

Air bubble at "3 o'clock" turn left foot clockwise, right foot counterclockwise.

Air bubble at "6 o'clock" turn both feet clockwise.


Air bubble at "9 o'clock" turn left foot counterclockwise, right foot clockwise.

- 2 The balance must be leveled and adjusted each time it is moved to a new location.




4.4 Assembling the Balance


- 1 Remove the transport protection (1).
 - 2 Insert the drip tray (2).
- ⇒ Push the tray in from the front over the bottom plate as far as the partition.



- Insert the grid weighing pan from the front.
- ⇒ Check that the grid weighing pan (1) (2) is correctly hooked in at both sides.




- 1 Insert the top door of the draft shield (1) at an angle (slightly less than 30 degrees) into the guide positioned at the **back**.
- 2 Swivel the draft-shield door (2) carefully down, **see** figure.




- The handles (A) must be turned toward the outside to allow installation of the side draft shield doors!


- 1 Insert the side doors of the draft shield according to the following instructions, **see** figure below.
- 2 Insert the side door at an angle of approx. 30° into the 2 openings, **see** figure.
- 3 Check that the side door is correctly inserted as shown!
- 4 Swivel the side door up against the balance until it engages with a click.
- 5 The side door must run easily, otherwise it is not correctly inserted.
- 6 Insert the second side door of the draft shield.
⇒ The procedure is identical.
- 7 Push the side doors completely to the back.




- 1 Insert the front glass (2) of the draft shield.
- 2 In the bottom part of the balance at the front, move at an angle from the top toward the bottom until the two hooks of the front glass of the draft shield lie on the rollers (1).
- 3 Swivel the front glass of the draft shield up until it engages.



- 1 Insert the terminal support.
- 2 First lay the cable in the guide by the terminal support.
- 3 Insert the terminal support into the opening of the front glass of the draft shield.
 - ⇒ The terminal support must engage with a click.



- 1 Mount the terminal.
- 2 Place the terminal in the center of the support.
- 3 Push it against the balance until it swivels slightly down at the front by the terminal support.
 - ⇒ You can push the cable into the balance.



Attention

The balance and the terminal are not fastened together by the terminal support! When transporting by hand, always hold the balance and the terminal firmly, **see** Transporting the Balance (page 17).

Note

You can also place the terminal free of the terminal support anywhere around the balance where the length of the cable allows.

4.5 Power Supply




WARNING

Risk of electric shock

- Ensure that the AC power pack for the balance is only used in accordance with the specifications listed in chapter General Data.
- Your instrument is supplied with a 3-pin power cable with an equipment grounding conductor. Only extension cables which meet this relevant standards and also have an equipment grounding conductor may be used. Intentional disconnection of the equipment grounding conductor is prohibited.


- Your balance is delivered complete with an AC adapter and a country-specific power supply cable. The AC adapter is suitable for all power supply voltages in the range of: 100 – 240 VAC, 50/60 Hz, for exact specifications, **see** Technical Data (page 24).
 - First, check whether the voltage of the power supply matches your local line voltage. If this is not the case, on no account connect the balance to the power, but contact your METTLER TOLEDO sales representative or dealer.
 - Guide the cables so that they cannot become damaged and will not be in your way during your daily work! Take care that the AC adapter cannot come into contact with liquids!
 - The power plug must be always accessible.
 - Before operating, check all cables for damage.
- Plug the AC adapter (1) into the socket (2) in the back of your balance and into the power supply.
- ⇒ After the balance has been connected to the power supply, it carries out a self test and is then ready for operation.



Note

If the display field remains dark, even though the power supply connection functions.

- 1 First disconnect the balance from the power supply.
- 2 Open the terminal.
- 3 Press both buttons (1) on the back of the terminal and open the upper part of the terminal.
- 4 Check that the plug for the terminal cable (2) is connected correctly inside the terminal.




4.6 Left/Right Operating of the Glass Draft Shield

The glass draft shield of your balance can be adapted to the environmental conditions and your weighing style, as well as to the type of weighing and loading.

The position of the handles determines which door(s) of the draft shield (left, right, or both) is/are opened.

Try various different combinations by moving the external handles into the upper or lower position. We recommend you to set up the glass draft shield so that it only opens on the side where the balance is loaded. Your balance then works faster, because there are fewer troublesome currents of air than when both doors of the draft shield are opened together.




4.7 Setting the Reading Angle and Positioning the Terminal

4.7.1 Optimise the Readability of the Terminal

Changing the reading angle

- 1 For a steeper reading angle, pull both levers (1) at the side upward.
⇒ You can then pull the upper part of the terminal slowly upward until it engages in the desired position. A total of 3 setting positions are available.
- 2 Move it in an appropriate position.
- 3 For a flatter reading angle, pull both levers (1) at the side upward, and press the lower part of the terminal downward.
- 4 Release both levers and the upper part of the terminal then engages in the desired position.



4.7.2 Remove Terminal and Place Close to the Balance


The terminal is connected to the balance by a cable. So you can arrange your workplace optimally, the terminal can be removed from the balance and placed separately.

Place the terminal separately

- 1 Switch the balance off.
- 2 Carefully lift the terminal off the terminal support.
You can leave the terminal support on the balance or remove it.
- 3 Pull the cable carefully out from the balance as far as possible.
- 4 Place the terminal where you want it to be.

Note

The cable can also be led out of the back of the balance. If working this way would be convenient for you, call your METTLER TOLEDO dealer who will adapt the balance for you.



4.8 Transporting the Balance

- 1 Switch off the balance.
- 2 The balance must be disconnected from the power supply.
- 3 Remove any interface cable from the balance.

4.8.1 Transporting Over Short Distances

If you wish to move your balance over a short distance to a new location, proceed as follows.




CAUTION

Damage of Device

Never lift the balance by the glass draft shield, as this can cause damage!

- 1 With one hand, hold the balance by the guide for the top door of the draft shield.
- 2 With your other hand, hold the terminal. The terminal is not rigidly fastened to the balance, so you must always hold the balance with one hand and the terminal with the other.
- 3 Carefully lift the balance and carry it to its new location, observe the notes in chapter Location (page 12).




4.8.2 Transporting Over Long Distances

If you want to transport or ship your balance over long distances, or if it is not certain that the balance will be transported upright, use the complete original packaging.

Disassemble the following parts

- 1 Lift the terminal (1) out of the terminal support and place it next to the support.
- 2 Pull the terminal support (2) off the balance.
- 3 Swivel the front glass (3) of the draft shield away from the balance.
- 4 Carefully fold the side doors (4+5) of the draft shield against the respective handles and pull the side doors out of the guide.
- 5 Swivel the front of the top door (6) of the draft shield up and pull the door out of the guide.
- 6 Carefully raise the front of the grid weighing pan (7) and lift it out of the guide.
- 7 Pull the drip tray (8) toward the front and out.




Pack the draft shield, the intermediate shelf and the terminal support (Pos. 2-6)

- Place these parts in the compartments provided in the original packaging.


Note

We advise you to place a sheet of paper between the sides glasses of the draft shield.



Pack the AC adapter, the power supply cable, and the individual parts (steps 7+8)

- 1 Place the AC adapter and the power supply cable in the packaging.
- 2 Place the drip tray (8) upside down in the packaging.
- 3 Place the grid weighing pan (7) upside down on the drip tray.
- 4 Place the ErgoClip "Basket" in the packaging.




CAUTION

Damage of Device

These instructions must be followed exactly, otherwise the balance may be damaged when inserting it into the packing cushions.


- 1 Push the transport protection over the weighing pan guide.
- 2 Push the guide of the top draft-shield door completely to the front.
- 3 Swivel the handle of the side door of the draft shield upward, and also push this door completely to the front.




Note

For packing both the balance and the terminal, you have a protective cover in which they were delivered. These are deliberately not shown in the illustrations so you can see better how the individual items must be positioned. However, we recommend you to use these protective covers.



- 1 Place the terminal on the balance (see illustration) and carefully insert the balance into the bottom packing cushion.
 - 2 Take the terminal and place it in front of the packing cushion on the table.
- Insert the packing set with the draft-shield glasses into the packing cushion, **see** illustration.




- Place the set with the AC adapter in front of the set with the draft-shield glasses.



- Insert the terminal into the packing cushion as shown in the illustration.




- 1 Now put the top packing cushion in place.
 - ⇒ Taking care to position it correctly.
- 2 Pass the lifting strap around both packing cushions, **see** illustration.
- 3 Tighten it until it lies close against the packaging.
 - ⇒ You can now lift the packed balance by the lifting strap and insert it into the transport carton.



4.9 Below-the-Balance Weighing

So that weighings can be carried out below the working surface (below-the-balance weighing), your balance is provided with a special hanger.

- 1 Switch off the balance.
- 2 Unplug the cable of the AC adapter from the back of the balance.
- 3 Also remove any interface cables.
- 4 Push all the doors of the glass draft shield completely to the back.
- 5 Lift the terminal off the terminal support.
- 6 Open the terminal and carefully pull the connecting cable out, **see** Power Supply (page 15).
- 7 Put the terminal down at the side of the balance.
- 8 Pull the balance over the table edge just far enough that you see the opening from below, **see** figure left.
- 9 Slacken the screw until the cover plate can be turned to the side and the hanger for weighing below the balance is easily accessible. You must now fasten the cover plate in the new position by tightening the screw, **see** figure right.
 - ⇒ Your balance is now ready for mounting your equipment for below-the-balance weighings.



4.10 Installing the ErgoClip

Attention

Before you install an ErgoClip you must switch off the balance «**On/Off**» key.

To install the ErgoClip included in the delivery, or an optional ErgoClip, please proceed as follows:

- 1 Remove the grid weighing pan (SmartGrid) from the balance.
- 2 Snap the ErgoClip onto the grid weighing pan.
- 3 Replace the grid weighing pan (SmartGrid) along with the installed ErgoClip.
 - ⇒ The optional "Flask" or "Tube" ErgoClips can be inserted directly, **see** Accessories (page 30).
- 4 Switch the balance on again «**On/Off**» key.



Important to know!

If you do not switch the balance off before you do the installation, the FACT function is not activated.

Reason

Addition of the ErgoClip causes the dead-load tolerance range of the balance to be exceeded. The balance therefore does not activate FACT, so as not to interrupt the **assumed** weighing operation.

When this status icon appears in the display, it means: "The balance wants to execute FACT" but cannot.



4.11 Installing the Single-use Aluminum Weighing Pan

Note

For standard operation with conventional tare containers, we do **not** recommend using this weighing pan. Its use may affect the stabilization time and degree of accuracy. The listed specifications are reached without a single-use weighing pan.




CAUTION

Hand injuries

Take care when handling the aluminum weighing pan, the corners and edges are extremely sharp!

- To install the single-use aluminum weighing pan, remove the grid weighing pan from the weighing chamber, **see** Assembling the Balance (page 13).

⇒ Only to be used for weighing highly-specialized tare containers.



- 1 Place the single-use aluminum weighing pan onto the grid weighing pan from above.
- 2 Fold the 4 side flaps under the bars of the grid weighing pan.

4.12 Installing the Grid Weighing Pan Cover

Note

For standard operation with conventional tare containers, we do **not** recommend using this weighing pan. Its use may affect the stabilization time and degree of accuracy. The listed specifications are reached without a weighing pan.



CAUTION


Hand injuries

Take care when handling the weighing pan, the corners and edges are extremely sharp!

Attention

With installed grid weighing pan cover, the balance does not switch to "Standby" mode!

- 1 For the installation, remove the grid weighing pan from the weighing chamber.
- 2 Gently press the cover onto the grid weighing pan.
- 3 Replace the grid weighing pan with the installed grid weighing pan cover.



5 Maintenance

5.1 Cleaning

Periodically clean the weighing pan, the drip tray, the housing, and the terminal of your balance using the brush supplied with it. The maintenance interval depends on your standard operating procedure (SOP).

Please observe the following notes



WARNING

Damage of balance

- The balance must be disconnected from the power supply.
 - Ensure that no liquid comes into contact with the balance, the terminal or the AC adapter.
 - Never open the balance, terminal or AC adapter – they contain no components, which can be cleaned, repaired or replaced by the user.
-



CAUTION

Damage of balance

On no account use cleaning agents which contain solvents or abrasive ingredients, as this can result in damage to the terminal overlay.

Cleaning

Your balance is made from high quality, resistant materials and can therefore be cleaned with a commercially available, mild cleaning agent.

- 1 To clean the weighing chamber thoroughly, swivel the glasses of the draft shield away from the balance and pull them out of their fastenings.
- 2 Carefully raise the front of the weighing pan and lift it out of the guide.
- 3 Pull the drip tray away from the balance.
- 4 When you replace these parts, make sure they are in the correct position.

Note

Please contact your METTLER TOLEDO dealer for details of the available service options. Regular servicing by an authorized service engineer ensures constant accuracy for years to come and prolongs the service life of your balance.

5.2 Disposal

In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

Thank you for your contribution to environmental protection.



6 Technical Data

6.1 General Data



CAUTION

Use only with a tested AC Adapter with SELV output current.
Ensure correct polarity

Power supply

Power supply connector with AC/DC adapter:

11107909

Primary: 100-240 VAC, -15%/+10%, 50/60 Hz

Secondary: 12 VDC \pm 3%, 2.0 A (with electronic overload protection)

Cable to AC adapter:

Design: 3-core, with country-specific plug

Note

Make sure the power supply plug is freely accessible

Power supply to the balance:

12 VDC \pm 3%, 2.0 A, maximum ripple: 80 mVDCpp

Protection and standards

Overvoltage category:

Class II

Degree of pollution:

2

Protection:

Protected against dust and water

Standards for safety and EMC:

See Declaration of Conformity

Range of application:

For use only in closed interior rooms

Environmental conditions

Height above mean sea level:

Up to 4000 m

Ambient temperature:

5-40 °C

Relative air humidity:

Max. 80% at 31 °C, linearly decreasing to 50% at 40 °C, non-condensing

Warm-up time:

At least **120** minutes after connecting the balance to the power supply; when switched on from standby-mode, the balance is ready for operation immediately

Materials

Housing:

Die-cast aluminum, lacquered, plastic and chrome steel

Terminal:

Die-cast zinc, lacquered, and plastic

Grid weighing pan:

Chrome-nickel steel X5CrNi18-10


6.2 Explanatory Notes for the METTLER TOLEDO AC Adapter

The certified external power supply which conforms to the requirements for Class II double insulated equipment is not provided with a protective earth connection but with a functional earth connection for EMC purposes. This earth connection IS NOT a safety feature. Further information about conformance of our products can be found in the brochure "Declaration of Conformity" which is coming with each product.

In case of testing with regard to the European Directive 2001/95/EC the power supply and the balance have to be handled as Class II double insulated equipment.

Consequently an earth bonding test is not required. Similarly it is not necessary to carry out an earth bonding test between the supply earth conductor and any exposed metalwork on the balance.

Because the balance are sensitive to static charges a leakage resistor, typically 10 k Ω , is connected between the earth connector and the power supply output terminals. The arrangement is shown in the equivalent circuit diagram. This resistor is not part of the electrical safety arrangement and does not require testing at regular intervals.



6.3 Model-specific Data

		XS64	XS104	XS204DR
Limit values				
Maximum capacity		61 g	120 g	220 g
Readability		0.1 mg	0.1 mg	1 mg
Tare range (from...to)		0 ... 61 g	0 ... 120 g	0 ... 220 g
Maximum capacity, fine range		–	–	81 g
Readability, fine range		–	–	0.1 mg
Repeatability (at nominal load)	sd	0.1 mg (60 g)	0.1 mg (100 g)	0.7 mg (200 g)
Repeatability (at low load)	sd	0.07 mg (10 g)	0.07 mg (10 g)	0.5 mg (10 g)
Repeatability, fine range (at low load)	sd	–	–	0.1 mg (10 g)
Linearity deviation		0.2 mg	0.2 mg	1 mg
Eccentricity deviation (test load)		0.15 mg (20 g)	0.3 mg (50 g)	0.3 mg (100 g)
Sensitivity offset (test weight)		0.9 mg (60 g)	1 mg (100 g)	1 mg (200 g)
Sensitivity temperature drift ¹⁾		0.00015%/°C	0.00015%/°C	0.00015%/°C
Sensitivity stability		0.0002%/a	0.0002%/a	0.0002%/a
Typical values				
Repeatability	sd	0.04 mg	0.04 mg	0.4 mg
Repeatability, fine range	sd	–	–	0.04 mg
Linearity deviation		0.1 mg	0.13 mg	0.3 mg
Eccentric deviation (test load)		0.06 mg (20 g)	0.15 mg (50 g)	0.16 mg (100 g)
Sensitivity offset (test load)		1.2 mg (60 g)	0.6 mg (100 g)	0.8 mg (200 g)
Minimum sample weight (according to USP)		120 mg	120 mg	1200 mg
Minimum sample weight (according to USP), fine range		–	–	120 mg
Minimum sample weight (U=1%, k=2)		8 mg	8 mg	80 mg
Minimum sample weight (U=1%, k=2), fine range		–	–	8 mg
Settling time		1.5 s	1.5 s	1.5 s
Settling time, fine range		–	–	1.5 s
Dimensions				
Balance dimensions (WxDxH)		263x453x322 mm	263x453x322 mm	263x453x322 mm
Weighing pan dimensions		78x73 mm (WxD)	78x73 mm (WxD)	78x73 mm (WxD)
Typical uncertainties and supplementary data				
Repeatability	sd	0.04mg +0.000015%·Rgr	0.04mg +0.00002%·Rgr	0.4mg +0.00005%·Rgr

		XS64	XS104	XS204DR
Repeatability, fine range	sd	–	–	0.04mg +0.00002%·Rgr
Differential linearity deviation	sd	$\sqrt{(40pg \cdot Rnt)}$	$\sqrt{(40pg \cdot Rnt)}$	$\sqrt{(120pg \cdot Rnt)}$
Differential eccentric load deviation	sd	0.00015%·Rnt	0.00015%·Rnt	0.00008%·Rnt
Sensitivity offset	sd	0.001%·Rnt	0.0003%·Rnt	0.0002%·Rnt
Minimum sample weight (according to USP)		120mg+0.045%·Rgr	120mg+0.06%·Rgr	1200mg+0.15%·Rgr
Minimum sample weight (according to USP), fine range		–	–	120mg+0.06%·Rgr
Minimum sample weight (U=1%, k=2)		8mg+0.003%·Rgr	8mg+0.004%·Rgr	80mg+0.01%·Rgr
Minimum sample weight (U=1%, k=2), fine range		–	–	8mg+0.004%·Rgr
Weighing time		4 s	4 s	3.5 s
Weighing time, fine range		–	–	4 s
Interface update rate		23/s	23/s	23/s
Usable height of draft shield		235 mm	235 mm	235 mm
Weight of balance		9.1 kg	9.1 kg	9.1 kg
Number of built-in reference weights		2	2	2
Weights for routine testing				
OIML CarePac		50 g F2, 2 g E2	100 g F2, 5 g E2	200 g F2, 10 g F1
	Weights	#11123003	#11123002	#11123001
ASTM CarePac		50 g 1, 2 g 1	100 g 1, 5 g 1	200 g 1, 10 g 1
	Weights	#11123103	#11123102	#11123101

sd = Standard deviation

Rnt = Net weight (sample weight)

Rgr = Gross weight

a = Year (annum)

1) In the temperature range 10...30 °C

		XS204	XS105DU	XS205DU
Limit values				
Maximum capacity		220 g	120 g	220 g
Readability		0.1 mg	0.1 mg	0.1 mg
Tare range (from...to)		0 ... 220 g	0 ... 120 g	0 ... 220 g
Maximum capacity, fine range		–	41 g	81 g
Readability, fine range		–	0.01 mg	0.01 mg
Repeatability (at nominal load)	sd	0.1 mg (200 g)	0.1 mg (100 g)	0.1 mg (200 g)
Repeatability (at low load)	sd	0.07 mg (10 g)	0.05 mg (10 g)	0.05 mg (10 g)
Repeatability, fine range (at low load)	sd	–	0.02 mg (10 g)	0.02 mg (10 g)
Linearity deviation		0.2 mg	0.2 mg	0.2 mg
Eccentricity deviation (test load)		0.3 mg (100 g)	0.3 mg (50 g)	0.3 mg (100 g)
Sensitivity offset (test weight)		1 mg (200 g)	0.8 mg (100 g)	0.8 mg (200 g)
Sensitivity temperature drift ¹⁾		0.00015%/°C	0.00015%/°C	0.00015%/°C
Sensitivity stability		0.0002%/a	0.0002%/a	0.0002%/a
Typical values				
Repeatability	sd	0.04 mg	0.04 mg	0.04 mg
Repeatability, fine range	sd	–	0.01 mg	0.01 mg
Linearity deviation		0.13 mg	0.13 mg	0.13 mg
Eccentric deviation (test load)		0.16 mg (100 g)	0.15 mg (50 g)	0.16 mg (100 g)
Sensitivity offset (test load)		0.8 mg (200 g)	0.4 mg (100 g)	0.6 mg (200 g)
Minimum sample weight (according to USP)		120 mg	120 mg	120 mg
Minimum sample weight (according to USP), fine range		–	30 mg	30 mg
Minimum sample weight (U=1%, k=2)		8 mg	8 mg	8 mg
Minimum sample weight (U=1%, k=2), fine range		–	2 mg	2 mg
Settling time		1.5 s	1.5 s	1.5 s
Settling time, fine range		–	3 s	3 s

		XS204	XS105DU	XS205DU
Dimensions				
Balance dimensions (WxDxH)		263x453x322 mm	263x453x322 mm	263x453x322 mm
Weighing pan dimensions		78x73 mm (WxD)	78x73 mm (WxD)	78x73 mm (WxD)
Typical uncertainties and supplementary data				
Repeatability	sd	0.04mg +0.000015%·Rgr	0.04mg +0.00002%·Rgr	0.04mg +0.00002%·Rgr
Repeatability, fine range	sd	–	0.01mg +0.00004%·Rgr	0.01mg +0.00003%·Rgr
Differential linearity deviation	sd	$\sqrt{(20\text{pg}\cdot\text{Rnt})}$	$\sqrt{(40\text{pg}\cdot\text{Rnt})}$	$\sqrt{(20\text{pg}\cdot\text{Rnt})}$
Differential eccentric load deviation	sd	0.00008%·Rnt	0.00015%·Rnt	0.00008%·Rnt
Sensitivity offset	sd	0.0002%·Rnt	0.0002%·Rnt	0.00015%·Rnt
Minimum sample weight (according to USP)		120mg+0.045%·Rgr	120mg+0.06%·Rgr	120mg+0.06%·Rgr
Minimum sample weight (according to USP), fine range		–	30mg+0.12%·Rgr	30mg+0.09%·Rgr
Minimum sample weight (U=1%, k=2)		8mg+0.003%·Rgr	8mg+0.004%·Rgr	8mg+0.004%·Rgr
Minimum sample weight (U=1%, k=2), fine range		–	2mg+0.008%·Rgr	2mg+0.006%·Rgr
Weighing time		4 s	4 s	4 s
Weighing time, fine range		–	6 s	6 s
Interface update rate		23/s	23/s	23/s
Usable height of draft shield		235 mm	235 mm	235 mm
Weight of balance		9.1 kg	9.1 kg	9.1 kg
Number of built-in reference weights		2	2	2
Weights for routine testing				
OIML CarePac		200 g F2, 10 g F1	100 g F2, 5 g E2	200 g F2, 10 g F2
	Weights	#11123001	#11123002	#11123001
ASTM CarePac		200 g 1, 10 g 1	100 g 1, 5 g 1	200 g 1, 10 g 1
	Weights	#11123101	#11123102	#11123101

sd = Standard deviation

Rnt = Net weight (sample weight)


Rgr = Gross weight

a = Year (annum)

1) In the temperature range 10...30 °C


6.4 Dimensions

Dimensions in mm.



6.5 Interfaces

6.5.1 Specifications of RS232C


Interface type:	Voltage interface according to EIA RS-232C/DIN 66020 (CCITT V24/V.28)	
Max. cable length:	15 m	
Signal level:	Outputs: +5 V ... +15 V (RL = 3 – 7 kΩ) –5 V ... –15 V (RL = 3 – 7 kΩ)	Inputs: +3 V ... 25 V –3 V ... 25 V
Connector:	Sub-D, 9-pole, female	
Operating mode:	Full duplex	
Transmission mode:	Bit-serial, asynchronous	
Transmission code:	ASCII	
Baud rates:	600, 1200, 2400, 4800, 9600, 19200, 384001) (firmware selectable)	
Bits/parity:	7-bit/even, 7-bit/odd, 7-bit/none, 8-bit/none (firmware selectable)	
Stop bits:	1 stop bit	
Handshake:	None, XON/XOFF, RTS/CTS (firmware selectable)	
End-of-line:	<CR><LF>, <CR>, <LF> (firmware selectable)	
	Pin 2: Balance transmit line (TxD) Pin 3: Balance receive line (RxD) Pin 5: Ground signal (GND) Pin 7: Clear to send (hardware handshake) (CTS) Pin 8: Request to send (hardware handshake) (RTS)	

6.5.2 Specifications of "Aux" Connection

You can connect the METTLER TOLEDO "ErgoSens" or an external switch to sockets "Aux 1" and "Aux 2". This allows you to start functions such as taring, zeroing, printing and others.

External connection







Connector:	3.5 mm stereo jack connector	
Electrical data:	Max. voltage	12 V
	Max. current	150 mA












7 Accessories and Spare Parts

7.1 Accessories

You can increase the functionality of your balance with accessories from the METTLER TOLEDO range. The following options are available:

	Description	Part No.
Printers		
	BT-P42 printer with Bluetooth connection to instrument	11132540
	Paper roll, set of 5 pcs	00072456
	Paper roll, self-adhesive, set of 3 pcs	11600388
	Ribbon cartridge, black, set of 2 pcs	00065975
	RS-P42 printer with RS232C connection to instrument	00229265
	Paper roll, set of 5 pcs	00072456
	Paper roll, self-adhesive, set of 3 pcs	11600388
	Ribbon cartridge, black, set of 2 pcs	00065975
	RS-P25 printer with RS232C connection to instrument	11124300
	Paper roll, set of 5 pcs	00072456
	Paper roll, self-adhesive, set of 3 pcs	11600388
	Ribbon cartridge, black, set of 2 pcs	00065975
	RS-P26 printer with RS232C connection to instrument (with date and time)	11124303
	Paper roll, set of 5 pcs	00072456
	Paper roll, self-adhesive, set of 3 pcs	11600388
	Ribbon cartridge, black, set of 2 pcs	00065975
	LC-P45 application printer with additional functions	00229119
	Paper roll, set of 5 pcs	00072456
	Paper roll, self-adhesive, set of 3 pcs	11600388
	Ribbon cartridge, black, set of 2 pcs	00065975
Optional interfaces		
	Second RS232C Interface	11132500


	Ethernet Interface for connection to an Ethernet network	11132515
	BT option: Bluetooth Interface for multipoint connection for up to 6 Bluetooth devices	11132530
	BTS option: Bluetooth Interface, single-point connection	11132535
	PS/2 option: Interface for connecting commercial keyboards and barcode readers	11132520
	LocalCAN option: Interface for connection of up to 5 LC (Local-CAN) instruments	11132505
	MiniMettler option: Interface MiniMettler, for connection to older (legacy) METTLER TOLEDO systems	11132510
	RS232 - USB converter cable – Cable with converter to connect a balance (RS232) to a USB port	64088427
Cables for RS232C interface		
	RS9 – RS9 (m/f): connection cable for PC, length = 1 m	11101051



RS9 – RS25 (m/f): connection cable for PC, length = 1 m


11101052

Cables for LocalCAN Interface



LC – RS9: Cable for connecting a PC with RS232C, 9-pin (f), length = 2 m

00229065



LC – RS25: Cable for connecting a printer or PC with RS232C, 25-pin (m/f), length = 2 m

00229050



LC – CL: Cable for connecting a device with METTLER TOLEDO CL interface (5-pin), length = 2 m

00229130



LC – LC2: Extension cable for LocalCAN, length = 2 m

00229115



LC – LC5: Extension cable for LocalCAN, length = 5 m


00229116



LC – LCT: Cable branch (T-connector) for LocalCAN

00229118


Cables for MiniMettler Interface



MM – RS9f: RS232C connection cable to MiniMettler interface, length = 1.5 m

00229029


Cables for Terminal



Terminal extension cable, length = 4.5 m

11600517


Cable, one-sided open (2-pin)



Cable between balance and AC adapter, length = 4 m


11132037

Auxiliary displays



BT-BLD Bluetooth auxiliary display for table mounting, 168 mm, LCD display with backlighting

11132555



LC/RS-BLD auxiliary display on bench stand, backlit (incl. RS cable and separate AC adapter)


00224200



RS/LC-BLDS auxiliary display for table or balance mounting, 480 mm, LCD display with backlighting

11132630

Sensors



ErgoSens, optical sensor for hands-free operation

11132601


LC-Switchbox



For connection of up to 3 balances with LocalCAN interface to a printer

00229220

Footswitches



Footswitch with selectable function for balances (Aux 1, Aux 2)

11106741



LC-FS foot switch with selectable function for balances with LocalCAN interface

00229060

Filling-Process Control



LV11 automatic feeder for automatic loading of small items on the balance

21900608

LV11 Draft shield door

11106715

SQC14 filling process control

Compact instrument with printer for control of up to 16 articles

00236210

Compact instrument with printer for control of up to 60 articles

00236211

Universal AntiStatic Kit



Universal antistatic kit complete U-shaped, including electrode and power supply

11107767

Optional: Second U-electrode*

11107764

Optional: Point-electrode*

11107765

* Power supply for optional, second U-elektrode

11107766

11107764, or for optional point-electrode

11107765

ErgoClips


Weighing kit for various weighing containers



ErgoClip Weighing kit

11106707

Delivery: 3 ErgoClips for Round-Bottom Flask, Weighing Boat and Tube, 20 Weighing Boat, 10 Single-use aluminium weighing pans.











ErgoClip "Basket"

11106747



ErgoClip "Titration Basket"

11106883

	ErgoClip "Weighing Boat"	11106748
	ErgoClip "Round-Bottom Flask"	11106746
	ErgoClip "small Flask"	11140180
	ErgoClip "Filter holder"	11140185
	ErgoClip "Stand"	11140170
	ErgoClip "Flask"	11106764
	ErgoClip "Tube"	11106784
	ErgoClip for Quantos	11141570



ErgoClip Syringe

30008288




ErgoClip Solution Kit

11140251



Single-use aluminium weighing pans, 10 units

11106711



SmardGrid Cover, chromium-nickel steel

11106709



Single-use weighing boats, 500 units

11106712



Grey drip tray


30038741



MinWeigh Door ideal for use with ErgoClip "Flask"


11106749

Density determination




Density kit

11106706



Sinker for density of liquids in conjunction with Density Kit
Calibrated (sinker + certificate)
Recalibrated (new certificate)

00210260
00210672
00210674



Calibrated thermometer with certificate

11132685

Pipette Calibration



Evaporation Trap, incl. adapter

11140043



Evaporation Trap large

11138440



1-channel suction pump complete
Hose 2 m for suction pump

11138268
11138132




Reagent reservoirs, 5 pcs.

11600616



Barometer


11600086




Precision thermometer with clip, not certified

00238767


Barcode Reader




RS232C Barcode Reader	21901297
The following accessories are needed for operation (not included):	
RS232 F cable	21901305
Null modem adapter	21900924
Plus one of the following:	
AC adapter 5 V for EU	21901370
AC adapter 5 V for US	21901372
AC adapter 5 V for GB	21901371
AC adapter 5 V for AU	21901370
	+ 71209966



RS232C Barcode Reader – Cordless	21901299
The following accessories are needed for operation (not included):	
Cradle	21901300
RS232 F cable	21901305
Null modem adapter	21900924
Plus one of the following:	
AC adapter 12 V for EU	21901373
AC adapter 12 V for US	21901375
AC adapter 12 V for GB	21901374
AC adapter 12 V for AU	21901373
	+ 71209966




PS/2 Barcode Reader, without cable	21901297
PS/2 wedge single cable	21901307



PS/2Y Barcode Reader, without cable	21901297
PS/2 wedge twin (Y) cable	21901308

Anti-theft devices



Steel cable	11600361
-------------	----------

Transport Cases



Transport case for analytical balances	11106869
--	----------

Protective Covers



Protective cover for terminal, "S" and "M" platform

11106870


Dust covers



Dust cover

30035838

Software



LabX Software for One Click™ Weighing Solutions

Enables you to perform One Click™ Standard Preparation, One Click™ Loss on Drying, One Click™ Sieve Analysis and many other applications.

Simply start the method with the One Click™ shortcut on the balance touchscreen. LabX guides you step-by-step through the SOP on the balance, performs your calculations automatically, and takes care of saving all your data. The complete solution can be tailored to match your process requirements. Visit www.mt.com/one-click-weighing for more information

Freeweigh.Net

on request


21900895

Various




Terminal and printer stand, mounting on balance

11106730



Wall fixture for terminal


11132665



AC/DC adapter (without power cable) 100–240 VAC,
50/60 Hz, 0.3 A, 12 VDC 2.25 A


11107909

Power cable CH	00087920
Power cable EU	00087925
Power cable US	00088668
Power cable IT	00087457
Power cable DK	00087452
Power cable GB	00089405
Power cable AU	00088751
Power cable ZA	00089728
Power cable BR	30015268
Power cable JL	00225297
Power cable IN	11600569
Power cable JP	11107881
Power cable TH, PE	11107880




IP54 protective housing for AC adapter

11132550



Level bubble mirror

11140150



Weighing table




11138042



Production Certificate "PRO"

11106895

7.2 Spare Parts

	Pos	Description	Part No.
	1	Side draft shield door	11106841
	2	Top draft shield door	11106842
	3	Front glass	11106843
	4	Grid weighing pan cover	11106709
	5	Grid weighing pan	11106333
	6	Foot screw	11106323
	7	Terminal support	11106539
	8	Drip tray	11106449
	Packaging complete	11106849	
	Export box	11106860	

8 Appendix

8.1 MT-SICS Interface Commands and Functions

Many of the instruments and balances used have to be capable of integration in a complex computer or data acquisition system.

To enable you to integrate balances in your system in a simple manner and utilize their capabilities to the full, most balance functions are also available as appropriate commands via the data interface.

All new METTLER TOLEDO balances launched on the market support the standardized command set "METTLER TOLEDO Standard Interface Command Set" (MT-SICS). The commands available depending on the functionality of the balance.

For further information please refer to the Reference Manual MT-SICS downloadable from the Internet under

► www.mt.com/xs-analytical

8.2 Procedure for Certified Balances

Preface

Certified balances are subject to the national, legal requirements of "non-automatic balances".

Switching on the balance

- **Switching on**
 - Immediately after being switched on, the balance displays 0.000.. g.
 - The balance is always started up with the "Factory setting" unit.
- **Switch-on range**
 - At maximum 20% of the type load, otherwise overload is displayed (OIML R76 4.5.1).
- **Stored value as switch-on zero point**
 - It is not permissible to use a stored value as a switch-on zero point; the MT-SICS M35 command is not available (OIML R76 T.5.2).

Display

- **Display of the weight value**
 - The "e" certification value is always shown in the display and is specified at the model designation plate (OIML R76 T.3.2.3 and 7.1.4).
 - If the display increment is lower than the "e" certification value, this is variably displayed for the net, gross and weighed tare. (Graying of the digits or certification brackets) (OIML R76 T.2.5.4 and 3.4.1).
- In accordance with guidelines, the tested display increment (certification value) is never lower than 1 mg (OIML R76 T.3.4.2).
- At balances with $d = 0.1$ mg, the digits below 1 mg are displayed in gray. These digits in brackets are printed. In accordance with legal metrology requirements, this illustration does not affect the accuracy of the weighing results.
- **Units of measurement**
 - The display and info unit are firmly set to g or mg (depending on the model).
 - The following applies for the "Custom unit":
 - No certification brackets.
 - The following names are blocked, this applies to upper and lower case letters.
 - All official units (g, kg, ct etc.).
 - c, ca, car, cm, crt, cart, kt, gr, gra, gram, grm, k, kilo, to, ton.
 - All names with "o" which can be replaced by a zero (Oz, Ozt etc.).

- **Identification of the weight display**

- Gross, net, tare and other weight values are accordingly marked (OIML R76 4.6.5).
 - Net for net when a tare value has been used.
 - B or G for gross.
 - T for the weighed tare.
 - PT for the specified tare.
 - * or diff for the difference between the net or gross.

- **Info field**

- The info weight value is handled metrologically in the same way as the weight value in the main display.

Printout (OIML R76 4.6.11)

- If a tare value is entered manually (PreTare), the PreTare value is always printed along with the net value (PT 123.45 g).
- The printed weight values are identified in the same way as the weight value on the display. I.e. N, B or G, T, PT, diff or *, with differentiation.

Example:

Single-range balance.

N	123.4[5] g
PT	10.00 g → for PreTare
G	133.4[5] g

DR balance with 100.00 g fine range.

N	80.4[0] g
T	22.5[6] g → for weighed tare
G	102.9[] g

Balance functions

- **Reset to zero**

- The zero range is limited to a maximum of $\pm 2\%$ of the full load (OIML R76 4.5.1).

- **Tare**

- No negative tare values are permitted.
- Tare immediate (TI) is not permitted, the MT-SICS `TI` command is not available (OIML R76 4.6.4).

- **1/xd**

- **e = d**

The 1/xd switchover is not permitted (OIML R76 3.1.2).

- **e = 10d**

This is only permitted in the case of the 1/10d switchover.

- **e = 100d**

Only the 1/10d and 1/100d switchover are permitted.

9 Index

A

AC adapter	24
AC Adapter	24
Accessories	30
Assembling the Balance	13
Aux Connections	29

B

Balance functions	43
Below-the-balance weighing	20

C

Certified balances	42
Changing the reading angle	16
Cleaning	23
Conventions	6

D

Dimensions	28
Display	42
Display field remains dark	16
Display of the weight value	42
Disposal	23

E

Environmental conditions	24
ErgoClip	20
ErgoSens	29

F

Features	5
----------	---

G

Glass draft shield	16
GLP	5
Good Laboratory Practice	5

I

Identification of the weight display	43
Info field	43
Installing the grid weighing pan cover	21
Interface	
MT-SICS	42
ISO 14001	5
ISO 9001	5

M

Materials	24
MT-SICS	42

O

Overview	9
----------	---

P

Pack	18
------	----

Position of the handles	16
Power supply	15, 24
Power supply voltages	15
Printout	43
Protection and standards	24

R

Remove terminal	16
Reset to zero	43
RS232C Interface	29

S

Safety information	7
Scope of delivery	11
Self test	15
Setting up	10
Single-use weighing pan	21
Spare Parts	41

T

Tare	43
Technical data	24
Transporting over long distances	18
Transporting over short distances	17
Transporting the balance	17

U

Units of measurement	42
Unpacking the balance	10

GWP® – Good Weighing Practice™

The global weighing guideline GWP® reduces risks associated with your weighing processes and helps to

- choose the appropriate balance
- reduce costs by optimizing testing procedures
- comply with the most common regulatory requirements

► www.mt.com/GWP

www.mt.com/excellence

For more information

Mettler-Toledo AG, Laboratory & Weighing Technologies

CH-8606 Greifensee, Switzerland

Tel. +41 (0)44 944 22 11

Fax +41 (0)44 944 30 60

www.mt.com

Subject to technical changes.

© Mettler-Toledo AG 05/2012

11781099A en

