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- ### ● How to Use the Object Marker
1. Attach the marker to the revolving nosepiece similarly as the objective.
 2. Remove the protective cap from the marker point.
 3. Bring the marker into the optical path that comes through the specimen area to be reobserved and press slightly the marker point on the specimen by means of the knurled ring. (Fig. 1)
 4. The inner diameter of the printed mark is 1.8mm . In advance of the reobservation of the specimen with the desired magnification, use an objective of $10\times$ or lower magnification to bring the printed mark into the view-field and relocate the object with ease.

Fig. 1
図1

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- ### ● How to Refill the Marking Ink
1. Take the marker off the revolving nosepiece.
 2. Hold the marker by the knurled ring near the marker point and unscrew the cap from the refill opening using a coin or a screwdriver. (Fig. 2)
 3. Refill the plastic inner tube with 2~3 drops of the stamp-er ink X for refilling (vermillion color) manufactured Shachihata Co.
 4. Keep the refilled marker upright for 12 hours or more to let the refilled ink reach the marker point.

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● Attention

1. This marker is to be used for the objective of the parfocal distance 45mm.
2. Refilled ink is durable for a long time, however, this merit may cause the mark wet and stain around the mark by touching with a finger, etc.
3. The mark can be clear away by wiping with alcohol.
4. Immersion oil, if any, on the specimen surface must be cleared away, otherwise the marking will be impossible.
5. Place the protective cap on the marker point for a long time storage.
6. Refill ink specified as Shachihata Stamper Ink X for refilling (Vermillion dye type) is available on order. Contact your dealer.

OBJECT MARKER Instructions

● How to Use the Object Marker

1. Attach the marker to the revolving nosepiece similarly as the objective.
2. Remove the protective cap from the marker point.
3. Bring the marker into the optical path that comes through the specimen area to be reobserved and press slightly the marker point on the specimen by means of the knurled ring. (Fig. 1)
4. The inner diameter of the printed mark is 1.8mm . In advance of the reobservation of the specimen with the desired magnification, use an objective of $10\times$ or lower magnification to bring the printed mark into the view-field and relocate the object with ease.

● How to Refill the Marking Ink

1. Take the marker off the revolving nosepiece.
2. Hold the marker by the knurled ring near the marker point and unscrew the cap from the refill opening using a coin or a screwdriver. (Fig. 2)
3. Refill the plastic inner tube with 2~3 drops of the stamp-er ink X for refilling (vermillion color) manufactured Shachihata Co.
4. Keep the refilled marker upright for 12 hours or more to let the refilled ink reach the marker point.

NIKON CORPORATION

Printed in Japan

36J/E 98.8.H.10