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Warnings

DO NOT FILL THE VAPORIZER WITH ANY ANAESTHETIC AGENT OTHER THAN THE ONE SPECIFIED ON THE FRONT LABEL. THE VAPORIZER IS DESIGNED FOR THAT AGENT ONLY. ANY OTHER AGENT THAN THAT SPECIFIED CAN PROVE TO BE DANGEROUS TO A PATIENT.

DO NOT CARRY THE VAPORIZER BY THE CONTROL DIAL.

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DO NOT IMMERSE THE VAPORIZER IN ANY LIQUID, INCLUDING WATER.

DO NOT STERILISE THE VAPORIZER.

DO NOT DRAIN THE ANAESTHETIC AGENT INTO ANY CONTAINER OTHER THAN A PROPERLY MARKED CONTAINER.

DO NOT FILL THE VAPORIZER UNLESS THE CONTROL DIAL IS IN THE 'OFF' POSITION.

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TO AVOID POTENTIAL LEAKS, BEFORE MOUNTING THE VAPORIZER ON THE SELECTATEC MANIFOLD, CHECK THAT THE MANIFOLD PORT VALVE 'O' RINGS ARE INTACT AND THAT THERE IS NO FOREIGN MATTER AROUND THE MATING SURFACES.

KEEP THE VAPORIZER UPRIGHT AT ALL TIMES.

DO NOT PUT WATER OR ANY OTHER SOLVENT IN A VAPORIZER

Cautions

Do not have the vaporizer serviced by anyone other than an Ohmeda Authorised Service Centre.

Turn the Vaporizer 'off' when not in use.

Before filling vaporizers with screw cap fillers, if the vaporizer is pressurised turn the screw cap slowly.

Before filling keyed filler vaporizers, if the unit is pressurised turn the screw slowly when removing the plug.

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1. Description

General

Tec 3 Vaporizers are designed for 'out of circuit' use in continuous flow techniques of inhalation anaesthesia.


The vaporizers are temperature and flow compensated so that their output remains relatively constant despite cooling due to vaporization and variations in inlet flow (see - 'Effects of Variables').

Each vaporizer is clearly labelled with the name of the anaesthetic agent for which it is designed.

Control Dial

A single control dial with a concentration scale calibrated in % of anaesthetic agent vapour per total volume (v/v) is employed to set the desired concentration of the anaesthetic agent.

To prevent accidental displacement of the control dial from the 'OFF' position to 'ON', a release button is incorporated in the dial assembly. Simultaneous depression of the release button and anti-clockwise rotation of the dial is necessary to achieve an 'ON' setting.

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2. Specification

Calibration

Vaporizers are calibrated at 21°C and at elevated temperatures. The variation in output with temperature, flowrate and duration of use is small, and the variation in output when used with Intermittent Positive Pressure Respiration is negligible.

The control dial is calibrated to the levels shown on the performance curves.

Resistance to Gas Flow

5 cm.wg at the 'OFF' setting at 5 litre/min O₂.

21-29 cm.wg. when delivering vapour at 5 litre/min O₂ at 22°C.

Duration of Use

The rate of consumption of anaesthetic agent depends primarily on flowrate and vapour output concentration. As an approximate working figure 1.0 ml of liquid anaesthetic is required to provide 200 ml of vapour.

The rate of evaporation of anaesthetic agent may with caution be used as an approximate method of checking that the delivered output is not grossly in error, and as a means of estimating how often the vaporizer is likely to need refilling.

The approximate hourly consumption of anaesthetic agents can be expressed as follows:

$$3 \times \% \times F$$

Where % represents the setting of the vaporizer output percentage, F represents the input flowrate in litre/min.

Example: If a vaporizer is set to deliver 2% at 6 litre/min total input gas flowrate:

$$\begin{aligned} \text{Approximate rate of agent consumption} \\ = 3 \times 2 \times 6 = 36 \text{ ml/hour.} \end{aligned}$$

The figures are intended only for clinical guidance and are approximate. They will vary depending upon the type of anaesthetic agent employed, accuracy of graduation of flowmeters etc. and will be grossly in error if the vaporizer drain port is not fully closed.

Liquid Capacity

Amount of anaesthetic agent to fully charge the vaporizer = 125 ml. (nominal)

Amount retained by Wick System = 35 ml. (nominal)

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

Fluotec 5%

EFFECT OF FLOWRATE 21-22°C Oxygen


EFFECT OF TEMPERATURE 5 litre/min Oxygen

Pentec


EFFECT OF FLOWRATE 23°C Oxygen

Flow = 5 litre/min

At the max setting the delivered concentration is the maximum percentage obtainable at the temperature of the vaporizer (i.e. the saturation concentration)


EFFECT OF TEMPERATURE 5 litre/min Oxygen

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








Tritec

EFFECT OF FLOWRATE
22°C Oxygen


EFFECT OF TEMPERATURE
5 litre/min Oxygen


Fluotec 8%

EFFECT OF FLOWRATE
22°C Oxygen


EFFECT OF TEMPERATURE
5 litre/min Oxygen

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








Enflurane 5%

**EFFECT OF FLOWRATE
22°C Oxygen**


**EFFECT OF TEMPERATURE
Flow 6 litre/min of 2:1 N₂O/O₂ mixture**


Enflurane 7%

**EFFECT OF FLOWRATE
22°C Oxygen**

**EFFECT OF TEMPERATURE
5 litre/min Oxygen**

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







Effects of Variables

Temperature

The effects of variation of temperature are normally negligible at commonly used combinations of dial setting and ambient temperature.

The vaporizer responds very slowly to change in ambient temperature and (to prevent the valve from closing completely) as a safety feature the temperature sensitive valve does not respond to temperatures below the range 12 - 15°C approx.

Should the vaporizer temperature be lower than this, then the output can be expected to be lower than that indicated on the dial.

At temperatures above the range shown on the performance curves, the vaporizer output may be unpredictably high - particularly if the temperature approaches the boiling point of the anaesthetic agent.

To avoid inaccuracies due to extreme temperatures the vaporizer should be allowed to attain a temperature in the range shown on the performance curves prior to use.

Pressure

Vaporizers are graduated in v/v percentage at 760 mm Hg. If the ambient pressure changes the v/v % will change so that at an ambient pressure P mm Hg the delivered percentage (D% v/v) -

$$\text{Equation 1} \quad D = \frac{\% \times 760}{P}$$

where % is the nominal setting of the vaporizer

It is generally accepted that the depth of anaesthesia depends on the inspired partial pressure of agent and not the concentration by volume of agent.

To obtain a consistent depth of anaesthesia when gross changes of barometric pressure occur it is necessary to change the v/v % in inverse proportion to the barometric pressure.

The vaporizer automatically does this and for practical clinical purposes the effects of the barometric pressure can be ignored.

Back Pressure Steady

Low and Moderate Pressures:

The vaporizer cannot distinguish between pressures at the outlet due to barometric pressure and pressures in excess of barometric due to steady back pressures applied by downstream components. Equation 1 therefore applies with term P now being the absolute pressure at the outlet (i.e barometric pressure plus back pressure). Steady back pressure reduces the v/v percentage.

High Pressures:

Pressures in excess of approximately 400 mm Hg should not be imposed on the vaporizer since these may overcome the loads imposed by internal thrust springs.

Low and Moderate Back Pressures:

Currently it is unlikely that the steady back pressure imposed by commonly used downstream components (other than some ventilators) will exceed 30 mm Hg at commonly used flowrates. Back pressures as high as this would reduce the delivered v/v percentage (at 760 mm Hg barometric pressure) to :-

$$\text{Equation 2} \quad \frac{760}{790} = 0.96$$

of what otherwise would be expected. Under normal clinical circumstances effects of this magnitude can be ignored.

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Some ventilators may impose higher steady back pressures of perhaps 100 mm Hg producing more significant depression of the v/v percentage. The increased patient uptake of agent with improved ventilation can often mitigate these effects so as to obviate the need to compensate for increased back pressure at the vaporizer.

High Back Pressures:

Pressures in excess of 400 mm Hg could conceivably occur during procedures similar to bronchoscopy or because of occlusion of downstream tubing and piping or for other reasons. These effects on v/v percentage cannot be precisely predicted but the most likely effects will be reductions in concentration (or small increases).

Back Pressure Fluctuating

Fluctuating back pressures may be imposed on the vaporizer by downstream components, and assisted or controlled ventilation to the patient. This can affect the vaporizer and increase the concentration by intermittently altering the pressures and hence the flow distribution within the vaporizer. The greatest effects are observed at combinations of very low flowrates and low dial setting with large and rapid pressure fluctuations and become progressively less important as the dial setting and flowrate increase and the magnitude and rate of cycling of the pressure fluctuations decrease.

In clinical use the vaporizers are considered unaffected by all fluctuating back pressures which would occur under all normal clinically encountered conditions appertaining to human anaesthesia.

Carrier Gas Composition

Small effects can occur when the carrier gas composition is changed from oxygen to air or nitrous oxide/oxygen mixture. As a general rule variation of output with carrier gas composition can be considered of negligible clinical significance since the effects, if any, are normally less than 10% of setting. Where changes occur the usual effect is that the output is slightly depressed when nitrous oxide is employed compared to the output when oxygen is the carrier gas.

The presence of nitrous oxide reduces the required inspired concentration of volatile agent and this mitigates this small depression in output from the vaporizer.

Time Out of Service

If the anaesthetic machine on which the vaporizer is fitted is left for a period of time with no gases flowing a concentration of agent may be observable at the machine outlet when the gas flow is turned 'on' and the vaporizers turned 'off'. This concentration can be expected to fall rapidly to zero (e.g. within about 15 seconds at 5 litre/min.). This phenomenon is a normal characteristic of anaesthetic vaporizers and anaesthetic machines.

Clinically this is normally considered to be insignificant because of the small volume of vapour involved.

Other Variables

Ambient temperature, input flowrate and duration often can affect delivered concentrations, particularly when the vaporizers are used at extremes of the usual clinical range.

The valve design and temperature compensation system of the Tec Vaporizers reduces the effects to levels when under most clinical conditions their effect on vaporizing performance is clinically not significant. The nominal performance characteristics should be consulted for further details.

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4. Principle of Operation

Vaporizer Sump and Valve Assembly

The duct system lies within the sump cover; above this the rotary valve, the lower face of which has ducts and a curved vapour control channel machined into its surfaces.

The spindle of the rotary valve is keyed to the control dial.

In the 'OFF' position, the rotary valve makes a direct link between the 'inlet' and 'outlet' within the vaporizer. When the dial is turned 'ON' the carrier gas is split into two streams, designated bypass and vaporizing chamber flow respectively.

The vaporizing chamber is lined by the two concentric wicks which enclose a nickel-plated copper helix so that the space is converted into a long spiral outlet channel. The wicks, dipping into the liquid, ensure that the vapour is maintained at saturation concentration in the gas leaving the vaporizing chamber.



The amount of agent picked up in the vaporizing stream will vary due either to variation in room temperature or to the cooling which takes place when the agent is vaporized. Each causes changes in the effective vapour pressures of the anaesthetic agents. Unless some form of compensating device was used, the output of the vaporizer for a given flow and dial setting would change with temperature.

The Tec 3 Vaporizer has a temperature compensating device (thermostat) which utilises a bi-metallic strip that deflects according to its temperature to control the proportion of carrier gas entering the vaporizing chamber.

If the temperature of the vaporizer falls, the thermostat closes and more carrier gas is allowed into the vaporizing chamber.

If the temperature of the vaporizer increases, the thermostat opens and less carrier gas is allowed into the vaporizing chamber.

In this way the output of the vaporizer remains constant under conditions of changing temperature.

Rotary Valve and Sump Assembly

1. Thermostat
2. Bypass Chamber
3. Wick Assembly
4. Rotary Valve
5. On
6. Vapour Control Channel
7. Vaporizing Chamber
8. Fresh Gas In
9. Direct link when Vaporizer is OFF
10. Fresh Gas Out
11. Flow Control (Vapour Channel)
12. Enriched Fresh Gas Out

Rotary Valve and Sump Assembly.
Drehventil- und Behälter-Baugruppe
Chambre d'évaporation et ensemble des valves

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5. Installation

The vaporizer may be incorporated into a back bar assembly, on a Selectatec Manifold, used on a table-top stand or mounted on a clamp fitted to the frame of the anaesthesia apparatus.

Back Bar

The vaporizer should always be mounted between the flowmetering unit and the patient breathing circuit — but upstream of any absorber or humidifier.

Check the integrity of the fittings to ensure that they are leak tight etc. If in doubt advice should be sought from the manufacturer of the equipment to which the vaporizer is attached.

Unless otherwise specified, all vaporizers are supplied as standard with 23 mm Cagemount fitting inlet and outlet tapers:

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

HANDLE WITH CARE

TO AVOID POTENTIAL LEAKS, BEFORE MOUNTING ON THE SELECTATEC MANIFOLD CHECK THAT THE MANIFOLD PORT VALVE 'O' RINGS ARE INTACT AND THAT THERE IS NO FOREIGN MATTER AROUND THE MATING SURFACES.

Mounting a Vaporizer

To Mount a Vaporizer

- 1) Check the mounting port valve 'O' rings are intact and not damaged, renew as necessary.
- 2) Make sure the vaporizer control knob is in the 'off' position.
- 3) Make sure the vaporizer locking lever is at the "unlock" position.
- 4) Carefully lower the vaporizer onto the manifold location so that the vaporizer's conversion block covers the two manifold port valves.
- 5) To lock the vaporizer onto the manifold, turn the locking lever to the "lock" position.

Mounting on a Selectatec Manifold

1. Unlock
2. Lock
3. Conversion Block
4. Check 'O' ring is intact before fitting vaporizer
5. Selectatec Manifold
6. Locking Lever

Mounting on a Selectatec Manifold
Montage am Selectatec Verteiler
Installation sur rampe Selectatec

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To Check for Proper Mounting:

- 1) If more than one vaporizer is fitted visually check to make sure the tops of vaporizers are level and at the same height. If a vaporizer is visibly out of line, remove it from the manifold and try to mount it again.
- 2) When the vaporizers appear to be level, are at the same height and have their locking levers in the locked position, attempt to lift each vaporizer from the manifold. If a vaporizer can be lifted off the manifold, it is not properly mounted. Try to mount the vaporizer again.
- 3) Check the machine for leaks with the vaporizer turned to 'O' (or lowest graduation) and at 'off'.

Vaporizer Removal

To Remove a Vaporizer:

- 1) Turn 'off' the vaporizer.
- 2) Turn the locking lever to the "unlock" position to free the vaporizer from the manifold.
- 3) Carefully lift the vaporizer straight up and off the manifold.

Base Ring

If a greater clearance is required below the vaporizer, the base ring can be removed. However, if the ring is removed, vaporizers fitted with Keyed Fillers will not stand upright on a shelf or bench and an appropriate means of mounting them upright should be provided such as separate base rings or a Selectatec Storage Rack.

The base ring is held in place by three grub screws located inside the ring.

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6. Operating Instructions

WARNING

KEEP THE VAPORIZER UPRIGHT AT ALL TIMES

DO NOT CARRY THE VAPORIZER BY THE CONTROL DIAL

HANDLE WITH CARE

Turning On

To turn the vaporizer 'ON' depress the control dial release button.

To avoid the inadvertent delivery of small concentrations the control dial should be turned to 'OFF' when the vaporizer is not in use.

Filling and Draining

General

WARNING

DO NOT FILL VAPORIZER WITH ANY AGENT OTHER THAN THE ONE SPECIFIED ON THE FRONT LABEL. THE VAPORIZER IS DESIGNED FOR THAT AGENT ONLY. ANY OTHER AGENT THAN THAT SPECIFIED CAN PROVE TO BE DANGEROUS TO A PATIENT.

DO NOT FILL VAPORIZER UNLESS THE CONTROL DIAL IS IN THE 'OFF' POSITION.

DO NOT TURN THE DIAL 'ON' DURING FILLING OR ATTEMPT TO FILL BEYOND THE 'FULL' MARK.


DO NOT DRAIN THE AGENT INTO ANY CONTAINER OTHER THAN A PROPERLY MARKED CONTAINER.

Periodically check the agent level. The vaporizer should be refilled at appropriate intervals. The vaporizer will function satisfactorily as long as there is agent visible in the sight glass.

The vaporizer should be filled and used in an upright position. Small deviations from the upright position will not affect the output or the safety of the vaporizer, but because the agent depth is shallow in relation to the diameter of the vaporizing chamber more frequent checks of the agent levels should be carried out when small deviations from the upright position occur in order to avoid obtaining a misleading impression of the amount of agent in the vaporizer.

It is good practice when the level is low, to drain the vaporizer of residual anaesthetic agent. The small amount of agent should be discarded.

At intervals - ideally not exceeding two weeks - the vaporizer should be drained into an appropriately marked container when the agent level is low and the agent discarded. Less frequent intervals may be used when the anaesthetic agent does not contain additives or stabilising agents.

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Screw Cap Filler

To Fill:

Dial should be in 'OFF' position. Remove filler cap by turning cap counter-clockwise. Be sure drain plug is closed.

CAUTION Turn cap slowly if vaporizer is pressurised.



Verify that the drug is the same as that specified on the front. Pour agent slowly into opening. Observe proper agent level through sight glass.

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Replace cap by turning cap clockwise. Cap should be tight to prevent leaks.

To Drain:

Remove the filler cap to reveal the drain plug. Using inverted filler cap unscrew the plug, but do not remove. Drain only into a properly marked container.

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

Keyed Filler

The filling system consists of three elements - the bottle collar, the adaptor and the filling and draining unit fitted to the vaporizer.

On later systems a simple ball valve is incorporated in the bottle adaptor to ensure rapid yet reliable filling, and should be used in preference to the earlier type.

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

To Fill:

- 1) Remove cap and seal from anaesthetic bottle. Check anaesthetic bottle neck is not chipped or damaged. Match Keyed Bottle Adaptor to Keyed Bottle Collar, and screw together until airtight. The bottle is then ready for filling Vaporizer.


- 2) Ensure that the vaporizer control is set in the 'OFF' position. Turn top retaining screw on Vaporizer Filler unit counter-clockwise and withdraw Filler Plug.


CAUTION Turn screw slowly if vaporizer is pressurised.


- 3) Grasp the Keyed Bottle Adaptor, and bottle, by the plastic tubing with the thumb on top of Adaptor at the keyway — with the two holes DOWNWARDS for filling. Fit into the Filler Socket normally sealed by the Filler Plug. ONLY the correct Adaptor will fit into this matching Filler Socket. Take care to bend tube slightly, so bottle is below inlet level to prevent spillage.


- 4) After insertion, turn the retaining screw clockwise and tighten to seal Filler Adaptor in the Filler Socket.

- 5) Raise bottle above level of Filler Socket, but avoid kinking tube. A steady stream of air bubbles should emerge from the inner tube within 2 seconds. If this does not occur perform procedures for clearing air locks.

- 6) When vaporizer is filled to the 'full' level, lower the bottle. Remove Bottle Adaptor. If any excess liquid drains from the Filler Socket, allow this to escape completely before inserting and tightening Filler Plug. (BE SURE to replace and tighten plug or gas will escape through Filler). The Vaporizer is then ready for use.

If the vaporizer is dry the level will fall slightly as the wicks absorb the agent.

To Fill

1. Keyed Bottle Adaptor
2. Keyed Bottle Collar
3. MUST BE AIRTIGHT
4. VAPORIZER CONTROL MUST BE AT 'OFF'
5. Retaining Screw
6. Filler Plug
7. Filler Socket
8. Liquid level indicates when full
9. LET EXCESS DRAIN BEFORE REPLACING FILLER PLUG

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

To Drain:

- 1) Insert the Bottle Adaptor in Drain (lower) Socket, with the two holes UP, and tighten retaining screw. ONLY the correct Adaptor will fit into the matching Drain Socket.


- 2) To let air vent from bottle and draining to occur, remove Filler Plug. Hold bottle slightly downward without kinking tube.


CAUTION Turn the screw slowly if vaporizer is pressurised.


- 3) Open Drain Valve on front counter-clockwise and Vaporizer will drain.

- 4) At completion of draining close Drain Valve (clockwise). If draining process stops see 'Procedures for clearing air locks'.

- 5) Remove the bottle and replace the Filler Plug in the Filler Socket. Tighten retaining screw (clockwise).

Procedure for Clearing Air Locks (Earlier Filling Systems)

Filling:

- Wait until air bubbles clear and slightly lower bottle (but keeping liquid level above filler socket) or
- Gently lower and raise bottle several times or
- Lower and remove bottle adaptor and repeat from step 3 or
- Slide conical liquid deflector on inner tube away from cap to cover side hole and repeat from step 3. Deflector must be repositioned correctly between cap and side hole immediately above side hole for draining.

Draining

Check conical liquid deflector is correctly positioned and reposition if necessary.

To Drain

- Drain Socket
- Retaining Screw
- Drain Valve

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Checking the Calibration

General

The clinical adequacy and reliability of Ohmeda Tec Vaporizers have been demonstrated by the most stringent proof of all - that of continued widespread satisfactory use throughout the world over many years involving hundreds of millions of successful administrations.

The performance of most Tec Vaporizers which are in clinical use is monitored by observing patient signs and consumption of anaesthetic agents. Some users may, however, wish to employ analysers to determine whether any abnormalities of performance have developed either as a routine procedure or as part of an investigation.

In order to achieve the reliability and consistency standards of Tec Vaporizers, Ohmeda use closely specified test conditions, test methods and detailed protocol in conjunction with training, experience and quality auditing systems. Because of this the full programme necessary to ensure that a vaporizer complies with Ohmeda standards cannot be practicably carried out in a field situation.

The following points should be considered when any measurements are being carried out on vaporizers to determine whether any abnormalities of performance have developed.

- 1) In order to predict the concentration the vaporizer can be expected to deliver the detailed nominal performance data, and the preceding comments should be taken into account.
- 2) The method of test used should not be such that it bears little relation to normal conditions of clinical use.
- 3) Any sampling techniques should be such as to ensure:

The sample is fully representative of the vaporizer output, which may not be a homogeneous mixture at the vaporizer outlet.

Absorption of agent by any connecting tubing is negligible.
- 4) If a number of vaporizers are being examined at the same time the probability of them all being consistently in error is so remote as to be negligible and the cause of any apparent error will lie in the test method employed.
- 5) Consistent and reproducible analytical techniques should be used.
- 6) If unexpected results are obtained it is a wise precaution to repeat the observations since the vaporizer may be more reliable than the techniques used to observe its performance.
- 7) If unexpected results occur it is also worthwhile checking for sources of error (e.g flowmeter, leaks, absorption by adjacent components etc.)
- 8) Full account should be taken of any extraneous effects on the analyser which may arise from changes in carrier gas composition.
- 9) If the anaesthetic machine on which the vaporizer is fitted is left for a period of time with no gases flowing, sensitive analysers may detect small concentrations of agent for a short time at the machine outlet after the gas flow is turned 'on' with the vaporizer turned 'off'. This concentration can be expected to fall rapidly to zero (e.g. within about 15 seconds at 5 litre/min.). This phenomenon is a normal characteristic of anaesthetic vaporizers and anaesthetic machines.

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- 10) At the O setting (where marked) it is not abnormal for small steady concentrations to be observable on sensitive analysers.

Analytical Techniques

For field checking of the state of calibration many techniques and analysers are available. Ohmeda would not recommend any one technique or analyser in preference to another but account must be taken of errors of use and calibration of analysers and their reliability must be realistically considered.

The following method of check can be used where special equipment is not available and where a secondary check on analysers is desired. The characteristics of the vaporizer are such that if the vaporizer is satisfactory at one dial setting it should be satisfactory at all other graduations.

- 1) Check that the vaporizer is full and has been at ambient temperature of 22°C for at least three hours.
- 2) With the vaporizer securely mounted open the drain until no more liquid will run out and close the drain.
- 3) Check that the dial is turned 'off' and carefully and quickly refill with a measured amount of agent without spilling. (Use about 70 ml). Close drain securely.
- 4) Leave vaporizer at nominally 22°C for one hour to ensure temperature has stabilised.
- 5) Set flowrate to 5 litre/min O₂.
- 6) Turn dial to 2%, note time and check that flowrate is still 5 litre/min - readjust as necessary.
- 7) Leave vaporizer at this setting for 30 minutes periodically check and adjust the flowrate as necessary. Turn vaporizer 'off' and turn 'off' O₂.
- 8) Drain as in 2) and measure amount liquid drained off.

The amount of liquid consumed in ml. should be as follows:

Fluotec	13.5
Enfluratec	15.5
Isotec	15.5

Appropriate measures to handle exhaust gases and spillage should be carried out.

It should be appreciated that the above check method is designed to be quick and easy under ordinary hospital conditions and that the method is somewhat imprecise. Nevertheless, it would be unusual for measured liquid consumption to vary from that given above by more than about 25%.

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7. Maintenance

WARNING

DO NOT MODIFY, TAMPER WITH OR DISASSEMBLE THE VAPORIZER BECAUSE OF THE DANGERS OF DAMAGING THE UNIT AND ALTERING THE ACCURACY OF GRADUATION.

Observation of the instructions given earlier, regular servicing and normal professional vigilance is normally all that is required to maintain the vaporizer in a safe working condition.

Schedule

Every Two Weeks:

The vaporizer should be drained into an appropriately marked container when the agent level is low and the agent discarded. Less frequent intervals may be used when the anaesthetic agent does not contain additives or stabilising agents.

Annually:

The vaporizer should be serviced at an Ohmeda Authorised Service Centre.

This service includes:

Complete disassembly of components.

Thorough cleaning.

Inspection for damage and wear.

Renewal of wicks, seals and damaged, worn or outdated items.

Lubrication where necessary.

Checks of the delivered vapour concentrations under closely defined conditions at different temperatures, and regraduation or adjustment where necessary.

Cleaning

WARNING

DO NOT PUT WATER OR ANY OTHER SOLVENT IN A VAPORIZER. A VAPORIZER SHOULD BE FILLED WITH THE SPECIFIED ANAESTHETIC AGENT ONLY.

Clean the exterior of the vaporizer with a damp cloth.

Never allow cleaning agents to accumulate in the filler, the gas inlet and outlet ports or around the control dial.

Contamination

If an incorrect volatile agent is put into the vaporizer:

- 1) Drain and discard all liquid.
- 2) Set the dial to 'Max' and scavenge the vaporizer with 5 litre/min O₂ until no trace of the contaminant can be detected.
- 3) Allow two hours for the vaporizer temperature to stabilize before proceeding with CAUTION.

If the contaminant is not volatile (e.g water), drain and return to Service Centre.

Repairs

Repairs should only be carried out by the manufactures' service representatives.

Note

It is not possible to change the graduation range or anaesthetic agent.

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The "OH" symbol and name Ohmeda, and the logo, is a service mark of the Ohmeda Group, Inc. and its subsidiaries.
Le symbole "OH" et le nom Ohmeda et le logo sont des marques commerciales de la Ohmeda Group, Inc. et de ses filiales.